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office@mediashop-group.com | www.mediashop.tv

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**EN****Product Contents:**

- 1x Gymform® SIX PACK EMS device
- 1x Gymform® SIX PACK for the abdominal area
- 1x Gymform® SIX PACK MINI for all muscles accessory
- 1x Gymform® belt SIX PACK
- 6x self-adhesive pads for the Gymform® SIX PACK for use on the abdominal area
- 2x self-adhesive pads Gymform® SIX PACK MINI for all muscles accessory
- 1x Nutritional guide
- 1x battery type CR2032

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**EN****Product description (Illustration 1):**

1. Gymform® SIX PACK EMS device
2. Button for switching on or increasing the intensity
3. Button to reduce the intensity or to switch off the display
4. Control light
5. Gymform® SIX PACK for the abdominal area
6. Gymform® SIX PACK MINI for all muscles accessory
7. Gymform belt® SIX PACK
8. Self-adhesive pads for the Gymform® SIX PACK for use on the abdominal area
9. Self-adhesive pads for the Gymform® SIX PACK MINI for all muscles accessory

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EN

## IMPORTANT NOTES

- Read the instructions before use. Keep the instructions in a safe place for future reference.
- Always give the instruction manual together with the product to a third party.
- Before starting a training program, you should consult your physician and physiotherapist. Gymform® SIX PACK is designed

for adults with good health. Nevertheless, the device must be used carefully. Gymform® SIX PACK moves the muscle tissue with electrical stimulation. It has pulse frequency controls and comes as close as possible to natural body and muscle activity with contraction and rest. No more than two treatments should be carried out per day, and no more than two days of rest between treatments. Maximum period of continuous use: 20 minutes.

- This equipment is not intended for use by children or by persons with impaired physical, sensory or mental abilities or by persons without appropriate experience or knowledge, unless they are supervised by a person responsible for their safety or receive the instructions necessary for the use of the equipment. Children must be supervised to prevent them from playing with the device. Do not allow children to clean and maintain the unit.
- The unit may cause personal injury if used improperly.
- This device is NOT a medical

device.

- This device is intended for domestic use only.
- The use of the device by persons belonging to one of the following groups is not recommended: elderly people, long-term patients, pregnant women, people with skin diseases, people with electronic implants (e. g. pacemakers), people with hypertension, people with cancer, epilepsy or cognitive dysfunction, people with fever, during breastfeeding, after a recent operation, during treatment for diabetes with insulin.
- The use of electrodes in the chest area can increase the risk of fibrillation of the heart muscle.
- Simultaneous use of a high-frequency electro-surgical device can cause burns at the contact points of the Gymform® SIX PACK and damage it.
- If a short-wave therapy device or a microwave is in operation in the immediate vicinity (e. g. 1 m), this can affect the performance of the Gymform® SIX PACK.
- Put away electronic devices such

as wristwatches, pedometers, mobile phones, etc. during use, as this may cause interference with the electronic device.

- Do not wear jewellery or accessories such as necklaces, belts or earrings when using the device, as this may cause incidents, injuries or damage.
- This device is a powerful electrostimulation device and can cause injuries if the training program is set to the highest level immediately. Always start at the weakest level and slowly increase the intensity level.
- Do not use when driving or sleeping.
- The device is not waterproof. Do not use in bathrooms or other rooms with high humidity.
- Although the device has been tested, the user is solely responsible for its use and consequences.
- Do not use the device on the neck, genitals, scars, moles, varicose veins or on the head. Do not apply to fresh scars, including inflamed skin, infected areas or areas where acne may develop; do not apply to

thrombosis or other areas with vascular problems (varicose veins) or other areas of the body with reduced sensitivity. Do not apply to areas with injuries or restricted mobility (fractures or sprains).

- Do not place the device directly on metal implants. Always apply to the nearest muscle.
- Do not connect buttons or pads to any other object.
- Do not make any modifications to the product.
- Only use electrodes that are included in the product contents.
- Do not glue 2 self-adhesive pads together.
- If one of the two pads has not been attached correctly, you will feel itching.
- Do not use defective, damaged or dirty pads as they could cause you burns or injuries.
- During use, make sure that the pads are in contact with the electrodes. If the electrode touches the skin directly, this can lead to injuries or burns.
- Clean the area on which the pad is to be used. The pad adhe-

res better to the skin when the hair has been removed in the appropriate area. Do not use the device on the same day as the hair has been removed.

- Do not use products such as body cream or lotion with the device, as this could result in poor adhesion.
- Apply the device to the desired training area and switch it on only after placing it on the body.
- Turn off the device before removing it from your body.
- If the skin or body does not feel normal in any way, stop using the device immediately and consult a doctor. You could otherwise suffer an accident or injury.
- If you use the device at its maximum capacity or for a long time, you can fatigue the muscles and cause prickling or reddening of the skin.
- Do not heat or cool the product before use as it may deform. Do not expose it to sudden changes in temperature during use, such as pouring hot water over it or using it in the sauna, etc.

**When using Gymform® SIX PACK, women must observe the following waiting times:**

- At least six months after birth (please consult a doctor beforehand).
- One month after inserting a contraceptive coil.
- At least three months after a caesarean section (please consult a doctor beforehand).
- The first days of menstruation.

**Inserting the battery:**

The Gymform® SIX PACK works with the supplied CR2032 battery pack, which is located on the back of the EMS device.

- Insert a coin into the groove in the cover and turn it to the left to open.
- Remove the protective foil from the battery. Insert the battery with the positive (+) pole facing up, turn the cover to the right to close it (see figure 2).

**ATTENTION:** As the battery power decreases, electrical stimulation becomes weaker. If the battery is flat, the indicator lamp will go out and the buttons may not work properly. Replace the battery with a new one. Before disposing the equipment, remove the batteries and dispose them safely. If the unit is not to be used for a long time, remove the batteries to protect it from damage.

**Use of Gymform® SIX PACK:**

1. Fasten the EMS device (1) to the Gymform® SIX PACK (5) using the push buttons. Make sure that the EMS device is properly attached.
2. Remove the plastic protection films from one side of the smaller self-adhesive pads (8) and attach them to the Gymform® SIX PACK. The pads must cover the grey area completely.
3. Remove the plastic protective film from the other side of the smaller self-adhesive pads (8) and place the product on the stomach.
4. Press the „+ / I“ button (2) once to turn the device on. The indicator lamp (4) will then flash and indicate that you are at intensity level 1.
5. Press the „+ / I“ button (2) again to gradually increase the intensity level. There are a total of 10 steps.
6. Press the „- / o“ button (3) to reduce the intensity levels. If you are at intensity level 1 and press the „- / o“ key, the device switches off.
7. We recommend starting with low intensity levels and gradually increasing the intensity as you familiarize yourself with the device.
8. To switch off the device immediately, press the „- / o“ button (3) for 2 seconds.
9. Gymform® SIX PACK has only one automatic program. The program switches off automatically after 20 minutes.
10. 10 minutes a day is enough. If you want even better results, use Gymform® SIX PACK as programmed 20 minutes daily. You can also combine Gymform® SIX PACK with physical activity (e. g. push-ups) for better success. To ensure that

the product is securely attached, you can use the Gymform® SIX PACK belt (7) and attach it to the Gymform® SIX PACK.

1. Place the Gymform® SIX PACK and turn on the unit. Select the desired intensity.
  2. Lay the belt (7) around the waist so that the buckle is on the outside.
  3. Fix the belt (7) so that it sits comfortably and the Gymform® SIX PACK rests firmly against the stomach.
11. Gymform® SIX PACK switches off automatically when the self-adhesive pads are not in contact with the skin or when the battery needs to be replaced.
12. After using the unit, reattach the plastic protective foils.

**IMPORTANT:** Never remove the 2 round orange-coloured labels on the back of the Gymform® SIX PACK. If one of these labels is not on the metal pins, replace them. If any of the 2 metal pins are exposed, do not use the product until the pins are covered again.

**Use of the Gymform® SIX PACK MINI:**

This element is designed for use on the following muscle groups: chest muscles, biceps and triceps, upper thigh front upper area, upper thigh back upper area, calves, buttocks, lumbar muscles, upper back area, shoulders

1. Place the EMS device (1) on the Gymform® SIX PACK MINI (6) with the aid of the push buttons.
2. Remove the plastic protection films from one side of the two larger self-adhesive pads (9) and attach them to the Gymform® SIX PACK MINI. The pads must cover the grey area completely.
3. Remove the plastic protective film from the other side of the self-adhesive pads (9) and place the product on the area to be treated (see figure 3). Replace the plastic protective film after using the device.

**NOTE:** WE RECOMMEND CHANGING THE PADS AFTER EVERY 30 APPLICATIONS. IF YOU NOTE THAT THE PADS ARE PREMATURELY NOT STICKING SO WELL, DAMPEN THE PADS WITH A SLIGHTLY WET SPONGE.

**Note:** When you start to feel heat at the place where the pads are placed, the pads are likely coming to the end of their operating life. Turn off Gymform® SIX PACK, remove the device from your body, attach new pads and restart the training session.

**Cleaning and maintenance**

DO NOT WASH Gymform® SIX PACK! WIPE IT WITH A DRY OR DAMP CLOTH.  
Wipe the self-adhesive pads with a finger moistened with water. Do not immerse in water! Store in a cool, dry place at room temperature. Do not peel off the insulating label.

**Troubleshooting**

**THE CONTROL LAMP DOES NOT LIGHT UP:**

- Cause:**
- The battery has not been inserted.
  - The battery has not been inserted correctly.
  - The battery is flat.

**Solution:**

- Insert the battery.
- Check that the positive (+) battery terminal symbol is visible on the upside of the battery.
- Change the battery.

**THE CONTROL LAMP LIGHTS UP BUT THE DEVICE DOES NOT OPERATE:**

- Cause:**
- No skin contact present.
  - The EMS device has not been properly attached to the Gymform® SIX PACK.
  - Intensity level is set at 0.
  - The microprocessor has deactivated itself.
  - The battery is almost empty.
  - The pads are not covering the grey area completely.

**Solution:**

- Make sure that the pad is in direct skin contact.
- Check that the console is properly attached to the Gymform® SIX PACK with the brackets.
- Press the + / I button to increase the intensity.
- Carry out a reset by removing the battery and replacing it after 1 minute.
- Change the battery.
- **Place the pads so that they cover the grey area completely.**

**THE BUTTONS DO NOT OPERATE CORRECTLY:**

- Cause:**
- The microprocessor has deactivated itself.
  - The battery is almost empty.
  - The pads are not covering the grey area completely.

**Solution:**

- Carry out a reset by removing the battery and replacing it after 1 minute.
- Change the battery.
- **Place the pads so that they cover the grey area completely.**

**DISPOSAL**

 Do not dispose this product at the end of its service life through normal domestic waste. It has to be disposed through a collection agency that recycles electric and electronic devices. This is pointed out by the symbol on the product, in the instruction manual, and on the packaging. Please inquire about local collection agencies that may be

operated by your distributor or your municipal administration. By recycling, making use of the materials, or other forms of recycling old units you are making an important contribution to protecting our environment.

Batteries are not part of your household waste and must be taken to an authorized collection point.

This product complies with the European directives.

#### TECHNICAL SPECIFICATIONS:

Initial waveform: Pulse in the form of a two-phase square wave  
Pulse duration: 320  $\mu$ s  
Pulse frequency: 3.3 ~ 33.3 Hz  
Matching voltage: 36 V +/-5 V load resistance 1000 Ohm  
Input voltage: 3 V  
Treatment duration: Automatic shut-off after 20 minutes  
Intensity: Adjustable from 1 to 10  
Mode: 1 variable frequency mode  
Operating conditions: 10°C~40°C (50°F~104°F); 30~85% RH  
Storage conditions: 10°C~50°C (14°F~122°F); 10~95% RH  
Maximum power consumption: 32mA

Power supply: 1 battery type 3V CR2032 (included in product delivery contents)

This symbol indicates that the device is protected against electrostatic discharge by a double insulation type BF.

Observe the operating instructions.

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